# MAS-100 NT®

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## USP 41 General Chapter <797>

### Pharmaceutical Compounding - Sterile preparations

#### Viable Air Sampling

- Evaluation of airborne microorganisms using volumetric collection methods in the controlled air environments shall be performed by properly trained individuals for all compounding risk levels.
- Impaction shall be the preferred method of volumetric air sampling.

#### Air Sampling Devices


- The instructions in the manufacturer's user manual for verification and use of electric air samplers that actively collect volumes of air for evaluation shall be followed.
- A sufficient volume of air (400–1000 liters) shall be tested at each location in order to maximize sensitivity.

#### Air Sampling Frequency and Process

- Air sampling shall be performed at least semiannually (i.e. every 6 months), as part of the re-certification of facilities and equipment for area where primary engineering controls are located.

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### Suggested Frequency of Sampling for Aseptic Processing Areas<sup>a</sup> (USP 41 <1116>)

Sampling Area/Location	Frequency of Sampling
<b>Clean Room/RABS</b>	
<i>Critical zone (ISO 5 or better)</i>	
Active air sampling	Each operational shift
Surface monitoring	At the end of the operation
<i>Aseptic area adjacent critical zone</i>	
All sampling	Each operating shift
<i>Other nonadjacent aseptic areas</i>	
All sampling	Once per day
<b>Isolators</b>	
<i>Critical zone (ISO 5 or better)</i>	
Active air sampling	Once per day
Surface monitoring	At the end of the campaign
<i>Nonaseptic areas surrounding the isolator</i>	
All sampling	Once per month

<sup>a</sup> All operators are aseptically gowned in these environments (with the exception of background environments for isolators). These recommendations do not apply to production areas for nonsterile products or other classified environments in which fully aseptic gowns are not donned.

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
## Spesifikasi teknis MAS-100 NT®

Informasi umum	
Flow rate	100 standard liters/min ± 2.5%
Standard sampling volumes	50, 100, 250, 500, 1000 liter
Freely definable sampling volumes	Volumes freely selectable between 1 and 2000 liters (volumes set to 0 are not displayed)
Perforated lid	
Material	Anodized aluminum
Dimensions	Ø 109 x90 mm
Weight	190 g
Lid geometry	300 holes, Ø 0.6 mm
Sampling efficiency	d <sub>50</sub> value 1.1 µm @ 100 SLPM (confirmed acc. ISO 14698)
Autoclaving	20 min at +121°C
Petri dish	For standard Petri dishes and contact plates

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Jenis Media	No. Katalog	Target Mikroorganisme
TSA settle plate - ICR triple bag	1.46001.0020	Total Bakteri
TSA settle plate with neutralizer		
LT - ICR triple bag	1.46050.0020	Total Bakteri
TSA settle plate with neutralizer		
LTHTh - ICR triple bag	1.46069.0020	Total Bakteri
TSA settle plate - LI single bag	1.46004.0020	Total Bakteri
TSA settle plate with neutralizer		
LTH - LI single bag	1.46015.0020	Total Bakteri
TSA settle plate with neutralizer		
LTHTh - LI single bag	1.46002.0020	Total Bakteri
SDA Settle plate - ICR Triple Bag	1.46577.0020	Total Jamur
SDA Settle plate with neutralizer		
LTHTh - ICR Triple Bag	1.46005.0020	Total Jamur
SDA Settle plate - LI Single Bag	1.46028.0020	Total Jamur
SDA Settle plate with neutralizer		
LTHTh- LI Single Bag	1.46052.0020	Total Jamur
SDA Settle plate with chloramphenicol - LI Single Bag	1.46003.0020	Total Jamur

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## Main unit:

Deskripsi	No. Katalog
MAS 100 NT® Air Sampler	1.09191.0001
MAS 100 NT Ex® Air Sampler	1.09194.0001

## Accesories :

Deskripsi	No. Katalog
MAS 100 NT® dust cover	1.09084.0001
MAS 100 NT® main charger (5 pin socket)	1.09200.0001
MAS 100 NT® main charger (DC Socket)	1.19133.0001
MAS 100 NT® extra perforated lid (300 x 0.6)	1.09195.0001
MAS 100® (old) extra perforated lid (400 x 0.7)	1.09088.0001
MAS 100 NT® tripod	1.09326.0001
MAS 100 NT® perforated lid for contact plates	1.09213.0001
MAS 100 NT® adapter for contact plates	1.09214.0001
MAS 100 NT® Battery pack Li-ion	1.09208.0001


## Referensi

Anderson, A. A. (1958)  
"New Sampler for the collection, sizing and enumeration of viable airborne particles.  
" Journal of Bacteriology, Vol 76, pp. 471-484.

Meier, R., Zingre, H. (2000)  
"Qualification of air sampler systems: The MAS-100®" Swiss Pharma 22, No. 1-2, pp. 15-21

Feller, W.  
"An Introduction to Probability Theory and its applications", John Wiley and Sons Inc.,  
New York, 1950

"Microbiological Control and Monitoring of Aseptic Processing Environments".  
U.S. Pharmacopoeia 41 , General Chapters <1116>, 2018

**PT Merck Chemicals and Life Sciences**  
Jl. TB Simatupang No. 8, Pasar Rebo  
Jakarta 13760, Indonesia  
[www.merckmillipore.com](http://www.merckmillipore.com) | [www.sigmal Aldrich.com](http://www.sigmal Aldrich.com) | [www.merck.co.id](http://www.merck.co.id)

Toll Free: 0800 140 1253 | (021) 2856 5600

SS\_Merck@merckgroup.com