

## Order Information

Tools & products for your Integrated safety right from the start

MERCK


Cat. No.	Product Description	Packaging
1.40000.9010	Extran® MA 05 liquid, alkaline, phosphate-free concentrate	10 L
1.40000.9025	Extran® MA 05 liquid, alkaline, phosphate-free concentrate	25 L
1.03834.0001	Fittings For Capillaries With 3.2 Mm O.D., For HPLC Bottle Adapter S40 (Pack of 10)	10 Units
1.03831.0001	HPLC bottle adapter S40 with 3 tube connections and 1 connection for exhaust air filter, solvents dis	1 Unit
1.03830.0001	HPLC-bottle adapter with 3 tube connections ID 3.2 mm, solvents supply	1 Unit
1.00801.0001	Label set acc. GHS, 80 x 110 mm	1 Set
Z108340-1PAK	Leather gloves, size L	1 pair
1.03835.0001	PTFE-ferrule for capillaries with 3.2 mm O.D., for HPLC bottle adapter S40 (pack of 10)	10 Units
9.20078.0001	Safety carrier for 2.5 L glass bottles	1 Unit
1.40140.0001	Safety carrier for 4 L glass bottles V2013	1 Unit
9.67213.0001	Safety Stand 2.5 L	1 Unit
B8276-200EA	Scienceware® biohazard disposal bags 25 in. (64 cm) x 35 in. (89 cm) x 0.002 in., red	200 pcs
B8151-100EA	Scienceware® biohazard disposal bags 8 1/2 in. (22 cm) x 11 in. (28 cm) x 0.00125 in., red,	100 pcs
Z108286-1PAK	Tempo-Shield heat-resistant gloves, Gloves, L 18 in.	1 pair
Z679267-50PR	Tyvek® disposable footwear, Overshoes, one size, white with elasticated top	50 Pairs
Z281409-1EA	Uvex® Stealth™ safety goggles gray frame	1 Unit

## My M Safety app - SDS & Label

Just download and use...

My M Safety app – your one-stop, mobile connection to Merck safety expertise around the world. The app allows users to find safety information on any product by snapping a photo of the 2D Barcode on the product label or searching for a substance. This includes:

- identifying a substance and any hazards it presents, to exposure control and first aid information
- download documents, such as CoAs and SDSs



My M Safety App



Cek harga?  
Yuk, tanya  
DR. M aja!



@merckindonesia



DR. M  
Asisten Virtual  
Merck Life Science



Scan me

**PT Merck Chemicals and Life Sciences**  
Jl. TB Simatupang No. 8, Pasar Rebo  
Jakarta 13760, Indonesia  
www.merckmillipore.com  
www.sigmaaldrich.com



Toll Free: 0800 140 1253



CS.Merck@merckgroup.com




@mercklifescience.id

# SAFETY



# INSIDE

Daily safety, risk prevention and emergency help that meet your high expectations.



The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

**Supelco®**  
Analytical Products

# SAFETY INSIDE

Trust our safety expertise



Scan me

We combine ingenious products and packaging solutions with individual support to keep you safe, simplify your lab work, and help you save resources. Discover daily safety, risk prevention and emergency solutions that match your high standards.

Discover more:  
[SigmaAldrich.com/lab-safety](https://SigmaAldrich.com/lab-safety)



## Safebreak bottles for acids

Our PE-coated glass bottles protect you from splinters and chemical burns in case of accidents



Integrated safety  
right from the start

## Risk prevention

Preventing harm starts with the right packaging and chemicals.



## Silica gel

Our orange gel contains an indicator based on iron-salt instead of the carcinogenic cobalt dichloride used in blue silica gel

## HF dosage bottle



This ingenious bottle combines the durability of HDPE with a secure, drip-free closure system for precise pouring and optimum protection against contact with hydrofluoric acid

## HDPE, metal & stainless steel containers

Safe handling of  
higher volumes



## HDPE bottles

Lightweight and unbreakable alternative to glass that offers easy handling



## Extran® detergents

Enjoy reliable, sustainable cleaning power with our nontoxic, biodegradable detergents

Tools &  
products  
for your



## Daily safety

Handling hazardous goods safely is always a challenge in daily lab work. We make it easy! With our perfectly matched components – chemicals, packaging, safety tools and withdrawal systems – you'll benefit from safe, convenient and contamination-free work even during highly sensitive applications.



## Withdrawal systems

Ingeniously safe solutions for:

- Solvents in bottles
- Solvents, acids & bases in higher volumes



## Swift, simple and efficient Emergency help

The fast, safe and easy way to clear up aggressive or unpleasant chemical spills



Check our  
**Chemizorb®**  
Special Promo!

## Bottle opening keys

Safe and easy handling from the moment you open a bottle



## Pouring aid

Our innovative aid simplifies pouring from 2.5 l glass bottles to protect your hands from chemical burns



## Safety carrier

The carrier features PE-foam bottle protection and a tightly closing screw lid to help you avoid contact with harmful chemicals, vapors and splinters



## Label set

Label your own mixtures with comprehensive safety information and GHS pictograms



## Chemizorb® absorbents

Order Information

Tools & products for your Integrated safety right from the start

Cat. No.	Product Description	Packaging
1.09996.0001	Adapter S 40	1 Unit
1.03832.0001	Air valve for HPLC bottle adapter S40	1 Unit
Z288497-1EA	Aldrich® Imperial mini-spill kits capacity 1.8 gal (6.8 L)	1 Unit
Z423270-1PAK	Best® N-DEX® laboratory gloves powder-free, thickness 4 mil (0.12 mm), size 8-9 L	100 pcs
B8401-1EA	Biohazard disposal stand Desktop size	1 pcs
Z401099-1EA	Biohazard spill kit Disposable, L x W x D 7 1/4 in. x 6 1/2 in. x 5 1/4 in	1 Unit
1.03836.0001	Blanking Plug For Capillary Connections With 3.2 Mm I.D., For Hplc Bottle Adapter S40 (Pack of 10)	10 Units
C2062-1EA	Bottle carrier, Nalgene® for 4 liter bottles	1 Unit
1.08801.0001	Bottle Key For Opening And Closing Bottles With S 40 And S	5 Units
1.01568.1000	Chemizorb® Granules Absorbent for spilled liquids	1 Kg
1.01568.5000	Chemizorb® Granules Absorbent for spilled liquids	5 Kg
1.01568.9021	Chemizorb® Granules Absorbent for spilled liquids	20 Kg
1.03874.2000	Chemizorb® H+ Absorbent and neutralizer for spilled acids, with indicator	2 Kg
1.03874.5000	Chemizorb® H+ Absorbent and neutralizer for spilled acids, with indicator	5 Kg
1.01591.1000	Chemizorb® HF Absorbent and neutralizer for spilled hydrofluoric acid, with indicator	1 Kg
1.12576.0001	Chemizorb® Hg Reagents and accessories for absorption of mercury	1 Set
1.01569.0001	Chemizorb® Hg Reagents refill pack for Cat. No. 1.12576.0001	1 Set
1.01596.1000	Chemizorb® OH- Absorbent and neutralizer for spilled alkalis, with indicator	1 Kg
1.02051.0500	Chemizorb® Powder Absorbent for spilled liquids	500 G
1.02051.9025	Chemizorb® Powder Absorbent for spilled liquids	25 Kg
Z408506-1PAK	Clavies® biohazard autoclave gloves gauntlet L 11 in. (28 cm), overall L 18 1/2 in. (47 cm)	1 pair
Z408492-1PAK	Clavies® biohazard autoclave gloves gauntlet L 5 in. (13 cm), overall L 13 in. (33 cm)	1 pair
Z183520-1PAK	Cryo-gloves® size L, water-resistant	1 pair
Z183512-1PAK	Cryo-gloves® size M, water-resistant	1 pair
Z677280-100EA	Dermatril™ nitrile gloves size L	100 pcs
Z677280-10X100EA	Dermatril™ nitrile gloves size L	10 x 100 pcs
Z677272-100EA	Dermatril™ nitrile gloves size M	100 pcs
Z677272-10X100EA	Dermatril™ nitrile gloves size M	10 x 100 pcs
Z677264-10X100EA	Dermatril™ nitrile gloves size S	10 x 100 pcs
Z677396-10X50EA	Dermatril™ P nitrile gloves size L	10 x 50 pcs

Order Information

Tools & products for your Integrated safety right from the start

Cat. No.	Product Description	Packaging
Z677396-50EA	Dermatril™ P nitrile gloves size L	50 pcs
Z677388-50EA	Dermatril™ P nitrile gloves size M	50 pcs
Z677361-50EA	Dermatril™ P nitrile gloves size S	50 pcs
Z187739-1KT	Dry-box glove/sleeve assembly, straight sleeves, neoprene gloves, port size 8 in., size 10	1 Kit
Z170062-1PAK	Economical neoprene dry-box gloves size 9, port size 8 in.	1 pair
1.03833.0001	Exhaust air filter HPLC bottle adapter	1 Unit
1.07558.9010	Extran® AP 11 powder, mildly alkaline	10 Kg
1.07558.2000	Extran® AP 11 powder, mildly alkaline	2 Kg
1.07558.9025	Extran® AP 11 powder, mildly alkaline	25 Kg
1.07563.9010	Extran® AP 12 powder, alkaline	10 Kg
1.07563.2000	Extran® AP 12 powder, alkaline	2 Kg
1.07563.9025	Extran® AP 12 powder, alkaline	25 Kg
1.07565.2000	Extran® AP 13 powder, alkaline with detergents	2 Kg
1.07565.9010	Extran® AP 13 powder, alkaline with detergents	10 Kg
1.07565.9025	Extran® AP 13 powder, alkaline with detergents	25 Kg
1.40006.2500	Extran® AP 17 liquid, alkaline concentrate	2.5 L
1.40006.9010	Extran® AP 17 liquid, alkaline concentrate	10 L
1.40006.9025	Extran® AP 17 liquid, alkaline concentrate	25 L
1.07559.2500	Extran® AP 21 liquid, acidic, concentrate (contains phosphoric acid)	2.5 L
1.07559.9010	Extran® AP 21 liquid, acidic, concentrate (contains phosphoric acid)	10 L
1.07559.9025	Extran® AP 21 liquid, acidic, concentrate (contains phosphoric acid)	25 L
1.07561.2500	Extran® AP 22 liquid, acidic concentrate (contains citric acid)	2.5 L
1.07561.9010	Extran® AP 22 liquid, acidic concentrate (contains citric acid)	10 L
1.07561.9025	Extran® AP 22 liquid, acidic concentrate (contains citric acid)	25 L
1.07570.2000	Extran® AP 41 Enzymatic 2 Kg	2 Kg
1.07570.9025	Extran® AP 41 Enzymatic 25 Kg	25 Kg
1.07555.9010	Extran® MA 01 liquid, alkaline, concentrate	10 L
1.07555.2500	Extran® MA 01 liquid, alkaline, concentrate	2.5 L
1.07555.9025	Extran® MA 01 liquid, alkaline, concentrate	25 L
1.07553.9010	Extran® MA 02 liquid, neutral, concentrate	10 L
1.07553.2500	Extran® MA 02 liquid, neutral, concentrate	2.5 L
1.07553.9025	Extran® MA 02 liquid, neutral, concentrate	25 L
1.40000.2500	Extran® MA 05 liquid, alkaline, phosphate-free concentrate	2.5 L