

# ZHE PACKAGED for Analyze Volatile Waste by TCLP

## Analyze Volatile Waste by TCLP

Do your waste samples contain volatile substances, such as isopropanol, acetone, or xylene? [Merck's zero head space extractor \(ZHE\) hazardous waste filtration system](#) has been designated by the USEPA as a suitable apparatus for TCLP.



Product Bundle	No. Cat.
ZHE Hazardous Waste Filtration System	YT30090HW
PRESSURE VESSEL 316SS 5L	XX6700P05
BALL VALVE 1/4"NPTF	YY2029348
VENT/RELIEF VALVE, 9 BAR	XX6700024
PRESSURE GAUGE SS 0-7BAR / 0-100PSI	YY1301015
QUICK RELEASE NIPPLE AND COUPLING, 1/4"	XX6700030
ZHE QUICK-DISCONNECT FEMALE SWAGELOK 1/PK	YT3009013
HEX NIPPLE, 1/4" NPTM 51MM, SS	XX6700125
TUBE ADAPTER SS 1/4" NPTM TO 1/4" TUBE	XX6702507
PP TUBING 6MM (0.25") OD x 3M	XX6702508
1/8 NPTF TO M-LUER NICKEL plated	XX3002567
ZHE TIGHT SYRINGE 1/PK	YT3009015
ZHE TEDLAR BAG 10/PK	YT3009016
ZHE TO TEDLAR BAG ADAPTER 1/PK	YT3009017
AP40 PREFILTER 90MM 100PK	AP4009000
ROTARY AGITATOR 115V/60HZ	YI320RAHW
PVC TUBING, 1/2" ID 10FT W/CLAMPS 125PSI	XX6700035
Forceps, blunt end, stainless steel 3pk	XX6200006P
CRM TCLP metal in soil	SQC005-225G
TCLP Metals in Soil – PT proficiency testing material	SPE005-225G

Discover more: [www.merckmillipore.com/tclp](http://www.merckmillipore.com/tclp)

## Contact Us



Contact	Mobile	Email
Hermanto	0815-8140-869	hermanto@merckgroup.com

PT Merck Chemicals and Life Sciences  
Arkadia Green Park - G building, 9th floor  
Jl. TB Simatupang No.Kav. 88, Pasar Minggu,  
Jakarta, 12520



CS.Merck@merckgroup.com



Scan to start  
Check Price &  
Quick Order

<https://bit.ly/Quick-Order-Merck>

# TCLP PACKAGED for Soil Sludge Analysis



## Toxicity Characteristic Leaching Procedure (TCLP)

Our products and services provide you a complete solution for successfully performing the TCLP method on your waste sample. TCLP determines the mobility of a specified list of organic and inorganic contaminants (metals, pesticides, herbicides, solvents) in liquid, solid, and multiphasic wastes. TCLP simulates leaching of these contaminants in a landfill environment over time. Before waste may be deposited in designated landfills, it must be tested by TCLP.

Product Bundle	No. Cat.
Hazardous waste filtration system 142 mm	YT30142HW
Gasket, TC 1 "1/2", OD 1 "1/2" Fluoroelast	YY2004057
PVC tubing, 1/2: ID 10 FTW/Clamps 125psi	XX6700035
Rotary agitator 115/60 HZ	YI320RAHW
AP40 142MM Disc 50/pk	AP4014250
GL Fibre Filter 142 MM 50/pk	APFF14250
Forceps, blunt end, stainless steel 3pk	XX6200006P
CRM TCLP metal in soil	SQC005-225G
TCLP Metals in Soil – PT proficiency testing material	SPE005-225G

Discover more: [www.merckmillipore.com/tclp](http://www.merckmillipore.com/tclp)

## Contact Us



Contact

Mobile

Email

Hermanto

0815-8140-869

hermanto@merckgroup.com

PT Merck Chemicals and Life Sciences  
Arkadia Green Park - G building, 9th floor  
Jl. TB Simatupang No.Kav. 88, Pasar Minggu,  
Jakarta, 12520



Scan to start  
Check Price &  
Quick Order

<https://bit.ly/Quick-Order-Merck>